
















The Right Thread for Every Suture

Our B. Braun Recommendations

Indication	Thread Structure	Article Number	Name	Sales Unit
 Soft Tissue Surgery, Parenchymatous Tissue	Monofilament, mid-term absorbable	C2022014	Monosyn® violet USP 4/0 (metric: 1.5), length: 70cm, needle: HR22, Racepack	Pack of 36
		C2022025	Monosyn® violet USP 3/0 (metric: 2), length: 70cm, needle: HR26, Racepack	Pack of 36
		C2022036	Monosyn® violet USP 2/0 (metric: 3), length: 70cm, needle: HR30, Racepack	Pack of 36
 Intracutaneous Suture	Monofilament, mid-term absorbable	C2022204	Monosyn® violet USP 4/0 (metric: 1.5), length: 70cm, needle: DS19, Racepack	Pack of 36
		C2022215	Monosyn® violet USP 3/0 (metric: 2), length: 70cm, needle: DS24, Racepack	Pack of 36
 Ligatures	Multifilament, mid-term absorbable	C0058209	Novosyn® violet 1 (metric: 4), length: 140cm	Pack of 36
 Abdominal Wall/ Linea Alba	Monofilament, long-term absorbable	C2024037	MonoPlus® violet USP 0 (metric: 3.5), length: 70cm, needle: HR30, Racepack	Pack of 36
		C2024058	MonoPlus® violet USP 1 (metric: 4), length: 70cm, needle: HR37S, Racepack	Pack of 36
 Skin Non-absorbable	Monofilament, non-absorbable	C0935204	Dafilon® blue USP 4/0 (metric: 1.5), length: 75cm, needle: DS19	Pack of 36
		C0935352	Dafilon® blue USP 3/0 (metric: 2), length: 75cm, needle: DS24	Pack of 36
 Callus	Multifilament, mid-term absorbable	G0068615	Novosyn® violet USP 8/0 (metric: 0.4), length: 15cm, needle: DLM6	Pack of 12
 Lashline	Monofilament, mid-term absorbable	C2022204	Monosyn® violet USP 4/0 (metric: 1.5), length: 70cm, needle: DS19	Pack of 36
 Hernias	Monofilament, non-absorbable	C3090041	Optilene® USP 3/0 (metric: 2), length: 75cm, needle: HR26, Racepack	Pack of 36
 Othematoma	Monofilament, non-absorbable	C0932655	Dafilon® blue USP 3/0 (metric: 2), length: 75cm, needle: GS60	Pack of 36
 Joint Capsule	Monofilament, long-term absorbable	C0024286	MonoPlus® violet 2/0 (metric: 3), length: 70cm, needle: HRT17 (M)	Pack of 36
 Dental	Monofilament, short-term absorbable	C0025007	Monosyn® Quick undyed 4/0 (metric: 1.5), length: 70cm, needle: DS19 (M)	Pack of 36
 Soft Tissue Surgery, Parenchymatous Tissue and Dental	Monofilament, mid-term absorbable	C2022004	Monosyn® violet USP 4/0 (metric: 1.5), length: 70cm, needle: HR17, Racepack	Pack of 36



Description				
Needle Curvature	Needle Shape	Needle Tip	Needle Length	Needle Size
S 1/8 circle	G Cutting, Golf Ball Technology	T Trocar Tip	Distance between reinforcement zone and needle tip, dimensions of long end in mm	f Thin-bodied Needle
V 1/4 circle	R Round Body	N Blunt Tip		ff Extra-thin-bodied Needle
D 3/8 circle	S Cutting	S Cutting Tip (Inside)		s Strong Needle Body
H 1/2 circle	L Lancet	C Short Cutting Tip		ss Very Strong Needle Body
F 5/8 circle		M Microtip		v Break-off Needle
G Straight		MP Micro Precision Tip		b Black Needle
P Asymptotic		m Micro Needle		TO Take-Off removable needle. With a straight tug of the needle holder, the needle can be removed from the suture fiber
J Hook		mV Microvascular Needle		
SK Ski Needle				

Example
HRT 37s = 1/2 circle round-bodied needle with trocar tip, length 37 mm, strong needle body

B. Braun Medical Ltd | VetCare | Thorncliffe Park | Sheffield | S35 2PW
Tel 0114 225 9000 | Fax 0114 225 9111 | www.bbraun-vetcare.co.uk

XX-SCR-05-22